

A Sample Assay Geometry for a Wide Range of Gamma Spectroscopy Sample Types and Volumes with a Single Efficiency Calibration and Still Achieve Reasonable Accuracy

2488749



Frazier Bronson CHP >> fbronson@canberra.com
Canberra Industries Inc., Meriden, CT 06450, USA

1 INTRODUCTION


- > Conventional "Laboratory Quality" gamma spectroscopy is capable of achieving very high quality accurate results
- > But there are many situations that do not demand this high quality
 - emergency response samples
 - samples from the initial and operational phases of decontamination projects
 - environmental remediation samples from the initial and operational phases
 - samples where regulatory compliance is not the primary purpose
 - samples which are expected to be well above or well below a decision value
- > For these situations, the important items frequently are:
 - getting the result quickly
 - minimal sample preparation time and labor
 - ability to easily handle a wide range of sample types
 - ability to easily handle wide range of sample sizes
 - minimal time spent preparing or choosing multiple efficiency calibrations

2 PROBLEM

- > With conventional "laboratory quality" gamma spectroscopy using radioactive sources
 - Calibrations usually done in advance
 - > Difficult to quickly adapt to new situations
 - Calibrations usually done for a few convenient matrices, e.g. water
 - > Estimations done to convert to proper density
 - Calibrations usually done for a few different sample geometries
 - > Labor and time must be spent to make the sample fit those geometries


4 COMPUTATIONAL EXPERIMENT DESIGN

- > Experiment performed to determine best sample-detector geometry
 - Sample chamber is 40cm diameter cylinder
 - Detector position choices are top or bottom - which is best?
 - Calibration choices are normal or massimetric - which is best?
 - ISOCS Uncertainty Estimator feature used to compute relative standard deviation with variable sample parameters
- Experiment 1
 - > Sample matrix was water
 - > Sample amount was random variable, all values from 5-40cm equally probable
 - > Both detector choices and both calibration methods tried
- Experiment 2
 - > Sample matrix was random variable, with all the following equally probable
 - Dry soil, cellulose, sand, concrete, mineralized soil, aluminum, plastic, 75%soil+25% iron
 - > Sample amount was random variable, all values from 5-40 cm equally probable
 - > Sample density was random variable, all values from 0.5-2.0 g/cc equally probable
 - > Both detector choices and both calibration methods tried




3 METHODS AND SOLUTIONS


- > Use mathematical calibrations, instead of radioactive sources
 - They are very quick to do and allow new situations to be easily accommodated
 - They work with any matrix and any density
 - > Concrete, steel, soil, air, vegetation, wood
 - ISOCS, a CANBERRA product, is widely accepted, and very versatile
 - Using these requires some skill by the operator for correct use



- Massimetric efficiency calibration
 - > Efficiency is the product of normal efficiency x mass of sample - i.e. counts per gamma per gram
 - > Once sample thickness is above a certain value, the efficiency is constant
 - > Example here shows normal and massimetric efficiency for bottle of water at energies from 20 to 1500 keV
 - Bottle on top of detector
 - Water filled from 10 to 30cm
 - > Massimetric efficiency almost constant for all fill heights
 - > Result automatically in activity/gram without weighing the sample



- > Use special geometries that are relatively invariant with sample type and volume
 - Count variable size samples with detector at a constant distance from opposite side of sample
 - > As sample size increases two competing effects
 - Bigger samples are closer to the detector, and higher efficiency
 - Bigger samples have more self-absorption, and lower efficiency
 - > Efficiency relatively constant with increasing sample size
 - > Same result for both spheres and cylinders
 - This geometry used in FastScan Whole Body Counters
 - One calibration for all sized people - see poster 2488744



5 RESULT AND CONCLUSION

- > Experiment 1 results - constant sample type, variable sample amount

- Top detector best with normal calibration
- Bottom detector best with massimetric calibration
- Bottom detector and massimetric calibration a little bit better, and since no weighing of the sample is required, is the preferred method.
- Standard deviation <15% due to calibration uncertainty

Det Loc'n	Fill Height	Fill Matrix	Density	Cal Type	20 keV	60 keV	200 keV	1000 keV	2500 keV
Top	5-40	Water	1	Normal	23	11	10	15	20
Bottom	5-40	Water	1	Mass	1	5	8	14	17

- > Experiment 2 results - variable sample type, variable density, and variable amount

- Top detector best with normal calibration
- Bottom detector best with massimetric calibration
- Bottom detector and massimetric calibration a little bit better, and since no weighing of the sample is required, is the preferred method.
- Standard deviation <30% due to calibration uncertainty

Det Loc'n	Fill Height	Fill Matrix	Density	Cal Type	20 keV	60 keV	200 keV	1000 keV	2500 keV
Top	5-40	Many	0.5-2.0	Normal	89	39	30	25	25
Bottom	5-40	Many	0.5-2.0	Mass	109	24	13	21	27

- > Conclusion

- Best sample geometry is with detector on the bottom
- Best calibration method is Massimetric calibrations
- Results with a standard deviation of <30% can be obtained with no sample preparation, and with a wide range of sample types and sample amounts



www.canberra.com